
HUMAN SOLUBLE FLT-1/VEGF-R1 ELISA KIT

FOR THE QUANTITATIVE DETERMINATION
OF HUMAN FLT-1/VEGF-R1
CONCENTRATIONS IN SERUM AND
PLASMA



FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

PURCHASE INFORMATION:

ELISA NAME	HUMAN SOLUBLE FLT-1/VEGF-R1 ELISA
Catalog No.	SK00114-06
Lot No.	
Formulation	96 T
Standard range	3.9-250 ng/mL
Sensitivity	100 pg/mL
Sample Volume	100 μΙ
Sample Type	Serum, Plasma
Dilution factor	Optimal dilutions should be determined by each laboratory for each application
Specificity	Human sFLT-1
Intra-assay Precision	4-6%
Inter-assay Precision	8-10%
Storage	2 °C-8 °C

Order Contact:
AVISCERA BIOSCIENCE INC.
2348 Walsh Ave., Suite C
Santa Clara, CA 95051

Tel: (408) 982 0300 Fax: (408) 982 0301

Email: Sales@AvisceraBioscience.com Info@AvisceraBioscience.com

www.AvisceraBioscience.com

INTRODUCTION

Human soluble FLT-1 immunoassay is a 3.5 - 4.5 hour solid phase ELISA designed to measure human FLT-1 in serum and EDTA plasma. It contains recombinant human FLT-1 and antibodies raised against this protein. It has been shown to accurately quantify recombinant human FLT-1. Results obtained with naturally occurring FLT-1 samples showed linear curves that were parallel to the standard curves obtained using the kit standards. These results indicate that the immunoassay kit can be used to determine relative mass values for natural human FLT-1.

PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for FLT-1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any FLT-1 present is bound by the immobilized antibody. After washing away any unbound substances, an antibody specific for FLT-1 is added to the wells. Following a wash to remove any unbound antibody reagent, an Anti-Rabbit IgG HRP conjugate is added to the wells. After washing away any unbound enzyme, a substrate solution is added to the wells and color develops in proportion to the amount of FLT-1 bound in the initial step. The color development is stopped and the intensity of the color is measured.

LIMITATIONS OF THE PROCEDURE

- _ FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.
- _ The kit should not be used beyond the expiration date on the kit label.
- _ Do not mix or substitute reagents with those from other lots or sources.
- _ It is important that the Dilution Buffer selected for the standard curve be consistent with the samples being assayed.
- _ If samples generate values higher than the highest standard, dilute the samples with Dilution Buffer and repeat the assay.
- _ Any variation in standard diluent, operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.
- _ This assay is designed to eliminate interference by soluble receptors, binding proteins, and other factors present in biological samples. Until all factors

have been tested in the immunoassay, the possibility of interference cannot be excluded.

MATERIALS PROVIDED

DESCRIPTION	CODE	QUANTITY
FLT-1 Microplate - 96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against FLT-1.	114-06-01	1 plate
FLT-1 Standard – 250 ng/vial of recombinant human FLT-1 in a buffered protein base with preservatives; lyophilized.	114-06-02	1 vial
Detection Antibody Concentrate – 105 μL/vial, 100-fold concentrated of an antibody against FLT-1 with preservatives; lyophilized.	114-06-03	1 vial
Positive Control - one vial of recombinant human FLT-1 in a buffered protein base with preservatives; lyophilized.	114-06-04	1 vial
Anti Rabbit IgG HRP Conjugate - 120 µL/vial, 100-fold concentrated solution of Streptavidin HRP conjugate	ARIGHRP	1 vial
Dilution Buffer - 12 mL of buffered protein based solution with preservatives	DB01	1 bottle
Antibody and HRP Diluent Solution - 30 mL of buffered protein based solution with preservatives	DB08	1 bottle
Wash Buffer - 50 mL of 10- fold concentrated buffered surfactant, with preservative.	WB01	1 bottle
TMB Substrate Solution- 11 mL of TMB substrate solution	TMB01	1 bottle
Stop Solution - 11 mL of 0.5M HCI	S-STOP	1 bottle
Plate Sealer	EAPS	1 piece
Plastic Pouch	P01	1 piece

STORAGE

Unopened Kit: Store at 2 - 8° C for up to 12 months. For longer storage, unopened Standard, Positive Control and Detection Antibody Concentrate should

be stored at -20 or -70 °C. Do not use kit past expiration date.

Opened / Reconstituted Reagents: Reconstituted Standard and Detection Antibody Concentrate Solution SHOULD BE STORED at -20 °C or -70°C for up to one month. ARIGHRP Conjugate 100-fold concentrated and other components may be stored at 2 - 8°C for up to 12 months.

Microplate Wells: Return unused wells to the plastic pouch with the desiccant pack. Microplate may be stored for up to 6 months at 2 - 8° C after opening.

OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Microplate shaker (250-300rpm).
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 100 mL and 500 mL graduated cylinders.

PRECAUTIONS FOR USE

All reagents should be considered as potentially hazardous. The stop solution contains diluted hydrochloric acid. Appropriate care, therefore, should be taken while handling this solution. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water.

SAMPLE COLLECTION AND STORAGE

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes before centrifugation for 15 minutes at $1000 \times g$. Remove serum and assay immediately or aliquot and store samples at \leq -20° C. Avoid repeated freeze-thaw cycles.

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤-20° C. Avoid repeated freeze-thaw cycles.

Note: Use Aprotinin (enzyme inhibitor) (Code No.: 00700-01-25) for ALL sample collection to prevent sample degradation. 0.5 TIU per ml of sample solution.

SAMPLE PREPARATION

Samples do not require dilution, but optimal dilutions should be determined by each laboratory for each application with a pretest.

Use polypropylene test tubes.


REAGENT PREPARATION

Bring all reagents to room temperature before use.

Wash Buffer - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 50 mL of Wash Buffer Concentrate into deionized or distilled water (450 mL) to prepare 500 mL of Wash Buffer.

FLT-1 Standard - Refer to vial label for reconstitution volume. Reconstitute the FLT-1 standard with 1 mL of Dilution Buffer. This reconstitution produces a stock solution of 250 ng/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Pipette 250 μL of Dilution Buffer into tubes #1 to #6. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 250 ng/mL standard serves as the high standard. The Dilution Buffer serves as the zero standard (0 pg/mL).

TUBE	STANDARD	DILUTION BUFFER	CONCENTRATION
Stock	powder	1000 μΙ	250 ng/ml
#1	250 μl of stock	250 µl	125 ng/ml
# 2	250 μl of 1	250 µl	62.5 ng/ml
#3	250 μl of 2	250 µl	31.25 ng/ml
#4	250 μl of 3	250 µl	15.6 ng/ml
# 5	250 μl of 4	250 µl	7.8 ng/ml
#6	250 μl of 5	250 µl	3.9 ng/ml



Detection Antibody - Reconstitute the Detection Antibody Concentrate with 105 μ L of **Antibody and HRP Diluent Solution (DB08)** to produce a 100-fold concentrated stock solution. Pipette 10.395 mL of Antibody and HRP Diluent Solution into a 15 mL centrifuge tube and transfer 105 μ L of 100-fold concentrated stock solution to prepare working solution.

ARIGHRP Conjugate - Transfer 120 μ L of 100-fold concentrated stock solution to 11.88 mL of **Antibody and HRP Diluent Solution** to prepare working solution. **Note:** 1x working solution of ARIGHRP Conjugate should be used within a few days.


Positive Control - Reconstitute the positive control with 1 mL of Dilution Buffer to make positive control solution.

ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that blank, standards, samples and positive control be assayed in duplicates.

- 1. Prepare all reagents and working standards as directed in the previous sections.
- 2. Remove excess micro-plate strips from the plate frame, return them to the plastic pouch with the desiccant pack.
- 3. Add 100 μL of Dilution Buffer to Blank wells (F4, F5).
- 4. Add 100 μL of Standard (from B2, B3 to G2, G3, and G4, G5), samples, or positive control (E4, E5) per well. Cover with sealer. Incubate for two hours on microplate shaker at room temperature. A plate layout is provided to record standards and samples assayed.

- 5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with 1x Wash Buffer (300 µL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
- 6. Add 100 μ L of Detection Antibody working solution to each well. Cover with sealer. Incubate for 2 hours on micro-plate shaker at room temperature.
- 7. Repeat the aspiration/wash as in step 5.
- 8. Add 100 μL of **ARIGHRP Conjugate working solution** to each well. Incubate for 60 minutes on micro-plate shaker at room temperature. **Protect from light.**
- 9. Repeat the aspiration/wash as in step 5.
- 10. Add 100 μ L of Substrate Solution to each well. Incubate for 8-12 minutes at room temperature on plate shaker. **Protect from light.**
- 11. Add 100 μ L of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green, or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
- Determine the optical density of each well within 15 minutes, using a micro-plate reader set to 450 nm.



CALCULATION OF RESULTS

Average the duplicate readings for each standard, positive control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a log-log curve fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw

a best fit curve through the points on the graph. The data may be linearized by plotting the log of the FLT-1 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

TYPICAL DATA

This standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

FLT-1 (NG/ML)	CORRECTED (450NM)
Blank	0 (0.128)
3.9	0.056
7.8	0.085
15.6	0.152
31.25	0.245
62.5	0.413
125	0.829
250	1.394

• Lot No.:

Positive Control: 19 – 40 ng/mL

CALIBRATION

This immunoassay is calibrated against a highly purified E Coli derived recombinant human VEGFRI/FLT-1(335-654).

SENSITIVITY

Twenty-five assays were evaluated and the minimum detectable dose (MDD) of FLT-1 was 100pg/mL.

SPECIFICITY

PROTEINS	CROSSREACTIVITY (%)
Human sFLT-1(335-654)	100
Human sFLT-1(32-327)	0
Mouse FLT-1/Fc	0
Human PIGF	0
Human sVEGF-RII	0
Human sVEGF-R III	0

SUMMARY OF ASSAY PROCEDURE

Add 100 μl of standard, samples, positive control to each well. Incubate for 2 hours on plate shaker at RT. Aspirate and wash 4 times. Add 100 μl Detection Antibody working solution to each well. Incubate for 2 hours on plate shaker at RT. Aspirate and wash 4 times. Add 100 μl ARIGHRP conjugate working solution to each well. Incubate 60 min on the plate shaker at RT. Protect from light. Aspirate and wash 4 times. Add 100 μl Substrate Solution to each well. Incubate 8-12 min on plate shaker. Protect from light. Add 100 μl Stop Solution to each well. Read 450nm

within 15 min